

Your value will be created by ALTOCHEM.

We have made a valuable history with our goods and minds. We hope to create more powerful history with you







A LTOCHEM is established in 2002 started as a chemical product manufacturing company supplying industrial adhesives including for medical use. Especially, not only manufacturing synthetic orthopedic casting tape & splint, but also we have been supplying all the raw materials, production lines and techniques to other manufacturers as a total service. With consistent and passionate R&D investment we have been keeping our product's quality superb, which brought us high reputations from all customers.

This made us set a firm footing as a top company in its field and we will lead a new way.



2014-2017

- Promising Export Company Award
- Inno-biz company Award
- Promising Small & Medium enterprise Award by Kyeonggi-Do

2012-2013

- Patent for High Heat Resistant Polymer Resin
- Goverment Procurement for the Korean Military Hospital

2010

- Patent for Casting Tape
- Patent for Cast Machine

2009

- FDA certified
- CE certified
- ISO 9001 certified
- ISO 13485 certified

2008

- Research institute certified
- Manufacture Medical Devices certified
- KGMP certified

2007

- Venture company certified
- 1,000,000 USD Export Award

2005

- Own facilities completion

2002



OBUSINESS AREA

Chemical / Electronic

- 01 | Water born Binder & Hardner
- 02 | Coating resin for special purpose gloves
- 03 | Electronic PSA (Pressure Sensitive Adhesive)
- 04 | Electronic heat emission material
- 05 | Vacuum casting resin(PU / Epoxy)
- 06 | Electronic binder
- 07 | Electronic Metal Paste & Ink

Medical

- 01 | Medical adhesives
- 02 | Synthetic orthopedic casting tape & splint
- 03 | Casting tape & Splint production machine
- 04 | Yarn(fiberglass)
- 05 | Fabric(fiberglass, polyester, non-woven)
- 06 | Resin & Additives
- 07 | Production techniques
- 08 | Dehumidification facilities
- 09 | ALTO CMC Thumb Brace™









ALTOCAST[™]

Synthetic orthopaedic casting tape

Historical Memory

20 years technical experiences / 25 Countries exporting recorde Korea military hospital supplier for 5 years in a row Overseas local site installation for manufacturing machine

• Fiberglass or Polyester Substrate

ALTOCAST™ is a synthetic orthopaedic casting tape made of knitted fiberglass or polyester substrate impregnated with water-activating polyurethane resin. The fabric is extensible in all directions, which makes it very comfortable in use.

Exposure to any moisture will cause the tape to become rigid. The setting time is approximately 4 minutes and it takes 20 minutes for full weight bearing.


Smooth Finish

High quality knitted substrate and special resin formula provide smooth surface and strong layer adhesion.

O Superb and Distingushed Color -

Excellent Moulding


Superior multi-direction stretch gives excellent moulding on curved area.







• Products Model Name



size	package	type	substrate	color
ALTOCAST™	•	•	•	• White Red
2in x 4yds (5.0cm x 3.6m)	•	•	•	Black Orange
3in x 4yds (7.5cm x 3.6m)	1 Roll/ Pouch	• Hand	•	Blue Yellow
4in x 4yds (10.0cm x 3.6m)	. 10 Rolls/ Inbox	Hard or Soft	Fiberglass or Polyester	Light Blue Violet
5in x 4yds (12.5cm x 3.6m)	. 100 Rolls/ Outbox	• 3011	Polyestel	. Green pink
iin x 4yds (15.0cm x 3.6m)	•	•	0	Light Green
	•	•	•	Graphic Shape
Altocast pouch packages	3 4 5	Altocast	c Color & packages	3

Adhesive Promoting Company

ALTO U-PADUnder Cast Padding



* Soft touch on skin / * Minimize allergic reaction due to no chemical treatment

Features

- 01 | Provides excellent water repellency and air circulation
- 02 | Elasticity and high thickness relieve the impact.
- 03 | Protection of the skin under casting tape from the growth of bacteria and odor
- **04** | Easy treatment There is no need to use a cotton bandage, so there is no hassle.
- 05 | Suitable for use with all casting tapes, synthetic, and plaster
- 06 | Available in a variety of widths / Applied to all parts of the body



ALTO U-PADUnder Cast Padding

Indications For Use













Precaution

- **1.** It should be applied according to the instructions of a doctor or specialist.
- 2. Do not use contaminated under padding
- 3. Store the product at room temperature with low humidity.
- 4. In case of problems, consult with a doctor or specialist

Storage

- 1. Store the product at room temperature with low humidity.
- 2. The shelf life is three years.

Product Sizes & Packaging

2 inch(Code no: AUP2100) / Size: 5.0cm x 1000cm (2 inch x 11 yds) / Packaging : 1 Roll/box **3 inch**(Code no: AUP3100) / Size: 7.5cm x 1000cm (3 inch x 11 yds) / Packaging : 1 Roll/box

4 inch(Code no: AUP4100) / Size: 10.0cm x 1000cm (4 inch x 11 yds) / Packaging: 1 Roll/box

5 inch(Code no: AUP5100) / Size: 12.5cm x 1000cm (5 inch x 11 yds) / Packaging : 1 Roll/box

6 inch(Code no: AUP6100) / Size: 15.0cm x 1000cm (6 inch x 11 yds) / Packaging : 1 Roll/box









ALTOSPLINT[™]

Synthetic orthopaedic Splint


Historical Memory

20 years technical experiences/ 25 Countries exporting record Korea military hospital supplier 5 years in a row Overseas local site installation for manufacturing machine

• Fiberglass/Non-woven substrate

ALTOSPLINT™ is a synthetic orthopaedic product intended for the emergency and temporary use. Its outside pads are non-woven fabric and the inside is 7-8 layers of fiberglass or a single layer of non-woven substrate impregnated with water-activating polyurethane resin.

Exposure to any moisture will cause the splint to become rigid. The setting time is approximately 4 minutes and it takes 20 minutes for full weight bearing.






Features

- 1. ALTOSPLINT does not require gloves for application.
- 2. ALTOSPLINT has superb x-ray permeability, which allows taking x-ray with the splint.
- 3. ALTOSPLINT makes it easier to shape the contours of the human body especially bony prominences such as joints, elbows, ankles and heels.
- 4. ALTOSPLINT's skin-contact-pad (unprinted side) is hydrophobic and non skin-contact-pad (printed side) is thin layer that allows to spray water on the printed side by a sprayer instead of dipping in water, which is easy to use and does not make patients feel unpleasant because of wet skincontact-pad.
- 5. ALTOSPLINT Roll type is economical since you can cut the roll as long as you require
- 6. ALTOSPLINT Pre-cut type is simple to use since it has various lengths.

• Free Needle-sprouting phenomenon

When cutting, fiberglass sprouts up and becomes painful needles but our unique heat treatment leaves no hazard.



size	package	type	substrate
2in x 5yds (5.0cm x 4.5m)	1 Roll / In box		
3in x 5yds (7.5cm x 4.5m)	1 Roll / In box	Name of the last o	
4in x 5yds (10.0cm x 4.5m)	1 Roll / In box	200 9 0	
5in x 5yds (12.5cm x 4.5m)	1 Roll / In box		
6in x 5yds (15.0cm x 4.5m)	6in x 5yds (15.0cm x 4.5m) 1 Roll / In box		
2in x 10in (5.0cm x 25cm)	10 EA / In box	Koll	
3in x 12in (7.5cm x 30cm)	10 EA / In box		Fiberglass or Non-woven
3in x 35in (7.5cm x 88cm)	5 EA / In box	1	fabric
4in x 15in (10.0cm x 38cm)	10 EA / In box	19 19 19 19 19 19 19 19 19 19 19 19 19 1	
4in x 30in (10.0cm x 75cm)	5 EA / In box	COLLEGE	
5in x 30in (12.5cm x 75cm)	5 EA / In box		
5in x 45in (12.5cm x 114cm)	5 EA / In box		
6in x 30in (15.0cm x 75cm)	5 EA / In box	One-step(Pre-cut)	
6in x 45in (15.0cm x 114cm)	5 EA / In box		



Adhesive Promoting Company

ALTO System

For synthetic orthopaedic casting tape & splint

• Historical Memory

20 Years technical experiences / 25 Countries exporting recorded Korea military hospital supplier for 5 years in a row Overseas local site installation for manufacturing machine

• The most advanced technologies

ALTO system for the production of casting Tape & Splint is operated by a touch panel, which is over 90% automated so less manpower is required and it is easy to control an uniform quality.





Basic information

Production capacity

casting tape machine - max 2,400 roll/day (8hrs) Splint machine - max. 1,000 m/day (10 hours)

Machine dimension (Width x Lenght x Height)

Casting tape machine – $3.5m \times 1.4m \times 1.8m$ Splint Machine – $18.0m \times 1.5m \times 2.7m$

Utility consumption

Electricity: 50 kw except air conditioning facilities

Air Pressure: 7-8 kgf/cm²

Required Site

Raw material storage: About 200m² Production Area: About 200m² Products Storage: About 100m²

* Additional spaces for office, laboratory, etc.

should be considered.



System components

Casting tape machine

Splint Machine

Air Compressor

Dehumidifying System

Raw material Drying/Heating System

Heating oven for resin

Polyurethane resin mixer & supplier

Sealing & printing system

Bobbin cutting system

Weighing instrument

Laboratory equipment

Required manpower

For the operation of casting tape machine

Machine operation and control -1 person

Pouch sealing - 1 person

Raw material handling and packing -1 person

For the operation of splint machine

Machine operation and control -1 person

Pouch sealing and control -1 person

Raw material handling - 2 person

Adhesive Promoting Company

ALTO Solution

For synthetic orthopaedic casting tape & splint

• Raw materials with the proven qualities

ALTOCHEM has been supplying the raw materials for the production of synthetic orthopaedic casting tape and splint to the majority of the manufacturers since 2002. Especially ALTO PU Resin and Fiberglass Tape are the most important key elements to the products quality that we have received high reputations from all users.



• Fiberglass/ Polyester Tape

ALTO fiberglass / polyester tape has excellent tensile strength and stretches very well in all directions that allows casting tape and splint high performance in moulding the contours of the human body especially bony prominences such as joints, ankles and heels.

ALTO fiberglass Tape has no needle-sprouting problem. When cutting, fiberglass sprouts up and becomes painful needles but our unique heat treatment leaves no hazard.



O Polyurethane Resin

ALTO PU Resin contains special oil instead of surfactant for lubrication. So it is not foamy and not sticky while applying, and its surface is very smooth after hardening.

ALTO PU Resin is designed to be hydrophilic that makes the water reaches deep into the resin. So it has super high compression strength and adhesion between the tape layers.







Technical Service

ALTOCHEM provides our client with full technical service such as Mass production, PU Resin mixing system, Dehumidification system, QC management, Maintenance.



• Etc-Raw materials & Part for ALTO Cast & Splint



Chemical & Electronic Products

ALTOCHEM has been manufacturing vacuum molding resin for over 10 years, focusing to develop Heat Emission film, PSA resin also. Finally we have succeeded in develop Heat Emission Film and Paste. We have completed developing Metal paste and adhesive for electronic part as new item and selling for electronic products.

Vacuum Casting Resin








3D Product

Silicone Mold

Vacuum Casting Products

• Raw Materials for Electronic & Special Resin


ALTOCHEM's electronic materials products have excellent electrical and mechanical properties. Our products are using for flexible PCB, LED, Power Module, EMI shielding.



Heat Emission



Special Adhesive & PSA Resin



Electronic & Special Resin

Other Products: ALTO CMC Thumb Brace™

The ALTO CMC Thumb Brace stabilizes the basal joint of the thumb (CMC-1) and places the thumb's metacarpal (metacarpal 1) in a functional position.



