

50g SMS Isolation Gown
(Non-Sterile & with Sealing Tape) Level3 & EN13795-1:2019

Products Specification

<ul style="list-style-type: none">◎Model No.: YD-JG◎Isolation Gown (Non-Sterile)◎Executive Standard: AAMI ANSI/AAMI PB70-2012 Level 2/3◎Material: 50g SMS Polypropylene fiber 100%◎Certificate: SGS Testing Report/FDA/CE◎Quality System: ISO 13485:2006◎Size: XL or Customized size◎Sealing: Ultrasonic with blue sealing tape	<ul style="list-style-type: none">◎Place of Origin:China◎Color: Blue◎Packing: 100pcs/case◎Supply Ability:150,000pcs per dayTreatment Available:<ul style="list-style-type: none">◎Anti-Static & Soft◎Absorbent / Hydrophilic◎Alcohol Repellent / Anti-bacterial
	



YADI ISOLATION GOWN
OPERATION INSTRUCTION

50g SMS Isolation Gown(Full Back)
Model No.: YD-JG
Size: XL
Level 3

【Intended Use】
The product provides general protection. It can be used as general isolation in clinic, ward or laboratory. It is not suitable to use in hospitals for infectious control purpose. It is not intended to use in operation room.

【Caution】
1. This product is for one-time use only. If the package is damaged, please do not use it.
2. It shall be properly treated as required and followed the local laws and regulations after use.
3. Never share your product with others.
4. Wash your hands thoroughly upon removal of the product.

【Operation instruction】
1. Pay attention to the package, production batch number and effective date before using this product.
2. Open the package and carry out routine operation according to the requirements.
3. handle according to the routine after use.

【Contraindication】 Use with caution when you are allergic to nonwoven fabrics.

【Storage conditions】
Packaged finished products should be stored in a warehouse with relative humidity less than 80%, conditions of storage: between -30℃ to +70℃, ventilated and dry, clean and hygienic, and free of harmful gases.

	End of shelf life: 2 years after manufacture date (Check the Approved Certificate for Product)
	Recommended conditions of storage: between -30℃ to +70℃
	Recommended conditions of storage: relative humidity <80%

☒ ☒ ☒

【Manufacturer】 Foshan Yadi clothing Co., Ltd.
【Add】 No. 4 building, yuanguzao porcelain factory, guzao village, Zhangcha street, Chancheng District, Foshan City China (528051)

YADI TEL: 86-757-83312999 / 86-757-83126811

Production Certificate
Level3

50g SMS Isolation Gown(Full Back)

MODELNO: YD-JG

EXECUTIVE STSNDARD: AAMIPB70

BRAND NAME: YADI

Size: (XL)

COMPONENT: 100% Polypropylene Fiber

QUANTITY: 100 PCS

INSPECTOR: QC01

LOTNO: 2020102001

PRODUCTION DATE: 2020-10-20

EXPIRY DATE: 2022-09-19

Foshan Yadi clothing Co., Ltd.
ADD: No. 4 building, yuanguzao porcelain factory, guzao village, Zhangcha street, Chancheng District, Foshan Guangdong, China 528051

The product has been inspected and meet the anticipated requirements, It is approved for shipment.

Art No.	Art Name	Component	Compliance With	Sewing Type	Color	Unit net/weight kg	Carton net weight kg	Carton Gross Weight KG	Carton Packing q'ty	Carton Dimension
YD-LT	Isolation Clothing (Non-Sterile) (One-Piece with Hood, without shoe cove)	45g SMS polypropylene fiber: 57.1% PE: 42.9%	Level 1 ANSI/AAMI PB70-2012	Ultrasonic	blue	0.2	20	21.3	100 pcs	60*42*37 CM
					white	0.15	15	16.3	100 pcs	60*42*37 CM
YD-G	Isolation Gown (non-sterile)	45 SMS polypropylene fiber: 55.5% PE: 44.5%	Level 2 ANSI/AAMI PB70-2012 EN13795-1:2019	Ultrasonic	blue	0.11	11	12.3	100 pcs	61*33*55CM
YD-G	Surgical Gown (non-sterile)	50g SMS Polypropylene fiber 100%	Level 3 ANSI/AAMI PB70-2012	Ultrasonic	blue	0.12	12	13.3	100 pcs	61*33*55CM
		45g SMS Polypropylene fiber 100%		Ultrasonic	blue	0.14	14	15.3	100 pcs	61*33*55CM
YD-JG	Isolation Gown (non-sterile) with blue sealing Tape	50g SMS Polypropylene fiber 100%	Level 3 ANSI/AAMI PB70-2012	Sewing+ blue Sealing Tape	blue	0.14	14	15.3	100 pcs	61*33*55CM

SGS Testing Report



中国认可
国际互认
检测
TESTING
CNAS L0599

Test Report

SL52045301313201TX

Date: November 06, 2020

Page 1 of 3

FOSHAN YADI CLOTHING CO., LTD
BLDG 4, GUZAO PORCELAIN, ZHANGCHA, ROAD, CHACHENG DISTRICT, FOSHAN CITY, GUANG DONG
PROVINCE, P.R.CHINA

The following sample(s) was/were submitted and identified on behalf of the client as:

Sample Description : (A) Isolation gown (non-sterile)

SGS Internal Ref.No. : GZF20-022982.001

Sample Color : (A) blue

Type No. : YD-JG

Manufacturer : FOSHAN YADI CLOTHING CO., LTD

Sample Receiving Date : Oct 28, 2020

Testing Period : Oct 28, 2020 - Nov 06, 2020

Test Result(s) : Unless otherwise stated the results shown in this test report refer only to the sample(s) tested, for further details, please refer to the following page(s).

Test Performed : Selected test(s) as requested by applicant

Conclusion: the provided sample meets ANSI/AAMI PB70-2012 Level 3.

Signed for and on behalf of
SGS-CSTC Standards Technical Services (Shanghai) Co., Ltd Testing Center

Sara Guo (Account Executive)



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are retained for 30 days only.
Attention: To check the authenticity of testing/inspection report & certificate, please contact us at telephone: (86-755) 8387 1443, or email: CN.Qcccheck@sgs.com
SGS-CSTC Standards Technical Services (Shanghai) Co., Ltd. 5th Building, No. 888, Yehen Road, Xuhui District Shanghai, China 200233 t (86-21) 61402666 f (86-21) 64958763 www.sgs.com.cn
中国·上海·徐汇区宜山路888号3号楼 邮编: 200233 t (86-21) 61402666 f (86-21) 64958763 e sgs.china@sgs.com

Member of the SGS Group (SGS SA)



Test Report

SL52045301313201TX

Date: November 06, 2020

Page 2 of 3

Test Result

Liquid Barrier Performance and Classification of Protective apparel and Drapes Intended for Use in Health Care Facilities

ANSI/AAMI PB70-2012

Section 4.2.1 Water Resistance: Impact Penetration Test
 (AATCC 42-2017)

As received

Weight of blotter gained (g)	1#	2#	3#	4#	5#
Area A (Critical zone-front)	0.0	0.0	0.0	0.0	0.0
Area B (Critical zone-sleeve)	0.0	0.0	0.0	0.0	0.0
Area C (Critical zone-Back)	0.0	0.0	0.0	0.0	0.0
Seam between areas A&B	0.0	0.0	0.0	0.0	0.0
Seam between areas A&C	0.0	0.0	0.0	0.0	0.0
Seam between areas B&C	0.0	0.0	0.0	0.0	0.0

Remark:

- Level 1: all critical zone components shall have a blotter weight gain of no more than 4.5grams(g)
- Level 2: all critical zone components shall have a blotter weight gain of no more than 1.0 grams(g)
- Level 3: all critical zone components shall have a blotter weight gain of no more than 1.0 grams(g)

Section 4.2.1 Water Resistance: Hydrostatic Pressure Test

(AATCC 127-2018; Hydrostatic Head; Rate of increase of water pressure: 60mbar/min; temp. of distilled water: 21°C, fabric face side of water)

As received

Water Column (cmH ₂ O)	1#	2#	3#	4#	5#
Area A (Critical zone-front)	>50.0	>50.0	>50.0	>50.0	>50.0
Area B (Critical zone-sleeve)	>50.0	>50.0	>50.0	>50.0	>50.0
Area C (Critical zone-Back)	>50.0	>50.0	>50.0	>50.0	>50.0
Seam between areas A&B	>50.0	>50.0	>50.0	>50.0	>50.0
Seam between areas A&C	>50.0	>50.0	>50.0	>50.0	>50.0
Seam between areas B&C	>50.0	>50.0	>50.0	>50.0	>50.0

Remark:

- Level 2: all critical zone components shall have a hydrostatic resistance of at least 20cmH₂O
- Level 3: all critical zone components shall have a hydrostatic resistance of at least 50cmH₂O



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are retained for 30 days only.

Attention: To check the authenticity of testing/inspection report & certificate, please contact us at telephone: (86-755) 8397 1443, or email: CN_Opscheck@sgs.com

3rd Building, No. 889, Yehon Road, Xuhui District Shanghai, China 200233 t (86-21) 61402666 f (86-21) 64958763 www.sgs.com.cn
 中国·上海·徐汇区宜山路889号3号楼 邮编: 200233 t (86-21) 61402666 f (86-21) 64958763 e sgs.china@sgs.com

Member of the SGS Group (SGS SA)



Test Report

SL52045301313201TX

Date: November 06, 2020

Page 3 of 3

Barrier performance of each component and final classification commended

	Impact Penetration Test AATCC 42(g)	Hydrostatic Pressure Test AATCC 127(cmH ₂ O)	Resistance to Bacteriophage Phi-X174 ASTM F 1671	Level	Final classification
Area A (Critical zone-front)	0.0	>50.0	/	Level 3	Level 3
Area B (Critical zone-sleeve)	0.0	>50.0	/	Level 3	
Area C (Critical zone-Back)	0.0	>50.0	/	Level 3	
Seam between areas A&B	0.0	>50.0	/	Level 3	
Seam between areas A&C	0.0	>50.0	/	Level 3	
Seam between areas B&C	0.0	>50.0	/	Level 3	

Remark:

- The barrier performance of all critical zone components, including seams and points of attachments, shall be determined. The classification of isolation gown shall be a number denoting the performance of the critical zone component having the lower barrier performance.
- Level 1: Impact Penetration Test-AATCC 42: ≤4.5g;
- Level 2: Impact Penetration Test-AATCC 42: ≤1.0g; Hydrostatic Pressure Test-AATCC 127: ≥20cmH₂O;
- Level 3: Impact Penetration Test-AATCC 42: ≤1.0g; Hydrostatic Pressure Test-AATCC 127: ≥50cmH₂O;
- Level 4: Resistance to Bacteriophage Phi-X174-ASTM F 1671: Pass.
- Only test the barrier performance of the sample. Label request, construction and other classed in ANSI/AAMI PB70-2012 is not checked.

Sample Photo



Face

Back

The statement of conformity in this test report is only based on measured values by the laboratory and does not take their uncertainties into consideration.

End of Report



SGS-CSR Technical Services (Shanghai) Co., Ltd.
Testing Center for Textiles

Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are retained for 30 days only.

Attention: To check the authenticity of testing/inspection report & certificate, please contact us at telephone: (86-755) 8387 1443, or email: CN.Doccheck@sgs.com


3rd Building No. 886, Yehon Road, Xuhui District Shanghai, China 200233 t (86-21) 61402666 f (86-21) 64958763 www.sgs.com.cn
中国·上海·徐汇区宜山路886号3号楼 邮编: 200233 t (86-21) 61402666 f (86-21) 64958763 e sgs.china@sgs.com

Member of the SGS Group (SGS SA)

Workshop Display



Ultrasonic Display





营 业 执 照

统一社会信用代码
91440604767347392Q

(副 本) (副本号: i-1)

 扫描二维码
或“国家企业信用
信息公示系统”了
解更多登记、备
案、许可、监管信
息。

名 称	佛山市雅迪制衣有限公司	注 册 资 本	伍拾万元人民币
类 型	有限责任公司(自然人投资或控股)	成 立 日 期	2004年10月19日
法 定 代 表 人	陶四平	经 营 期 限	长期
经 营 范 围	加工、制造：服装、床上用品；销售：服装、纺织原料（不含危险化学品），纺织品、羽绒、羽绒制品；研发、销售：熔喷布；研发、制造、销售：第一类、第二类、第三类医疗器械，卫生材料及医药用品，无纺布。（依法须经批准的项目，经相关部门批准后方可开展经营活动。）		
住 所	佛山市禅城区张槎街道古灶村原古灶瓷厂内4号楼（住所申报）		

登 记 机 关

2020 年 06 月 28 日



国家企业信用信息公示系统网址: <http://www.gsxt.gov.cn>

国家市场监督管理总局监制

FDA

 **U.S. FOOD & DRUG**
ADMINISTRATION

Follow FDA | En Español

SEARCH

Home Food Drugs Medical Devices Radiation-Emitting Products Vaccines, Blood & Biologics Animal & Veterinary Cosmetics Tobacco P

Establishment Registration & Device Listing

● FDA Home ● Medical Devices ● Databases

1 result found for **Owner Operator Number :**
10076526

New Search

Establishment Name	Registration Number	Current Registration Yr
FOSHAN YADI CLOTHING CO., LTD CHINA	3017417657	2021
Non-Surgical Isolation Gown - Surgical Gown(Non-Sterile) YD-LY, YD-Lg		Manufacturer
Suit, Surgical - Isolation Clothing(Non-Sterile) YD-LT		Manufacturer



QCC ISO Certificate No.: QCC638ISOQ



European & USA Size Comparing Chart

Size	Measuring (unit: inch)											
	A	B	C	D	E	F	G	H	M	N	J	K
UL	27	35	47	10	2 1/4	43	70	7	1 3/8	2 1/2	2 1/2	1 1/2

XL	28	37	48	10 1/2	2 1/4	44	70	7	1 3/8	2 1/2	2 1/2	1 1/2
Tolerance	±1	±1	±1	±1	±1/4	±1	NA	±1	±1/2	±1/2	±1/8	±1/4


Guangzhou King Yon Medical Tech Co.,Ltd

Address: No.718, Blv San Yuan Li, Bai Yuan
District, Guang Zhou City, Guang Dong
Province, P.R.China

Post Code: 510403

Website: www.King-Yon.com (Head Office)

Attn: Sam Deng (Mr.)

G.Manager

E-Mail: sales@King-Yon.com

Mobile Phone: +86-15089203268

WeChat ID:


Foshan Yadi Clothing Co.,Ltd

Address:Bldg 4,Guzhao Porcelain, Zhangcha
Road, Chancheng District, Foshan City,
Guang Dong Province, P.R. China

Post Code: 528051

Website: www.YadiSuits.com (Manufacture base)

Attn: Zeego (Mr.)

Sales Manager

E-Mail: sales@King-Yon.com

Mobile Phone: +86-13802807828

WeChat ID:

