

L.ma LED SKIN CARE MASK LED PRO VITSRO



GYEONGNAM HEADQUARTERS
71-107, Sugok-ro, Naeseo-eup, Masanhoewon-gu, Changwon-si, Gyeongsangnam-do, Republic of Korea / T : 055-238-3331

SEOUL BRANCH The 2nd floor, 14, Neungdong-ro 3-gil, Gwangjin-gu, Seoul, Republic of Korea / T : 02-469-3331

CHUNGCHEONG BRANCH 39-51, Buldang 4-ro, Seobuk-gu, Cheonan-si, Chungcheongnam-do, Republic of Korea / T : 041-556-6626

DAEGU BRANCH 21, Daehak-ro 32-gil, Gyeongsan-si, Gyeongsangbuk-do, Republic of Korea / T : 053-815-6627

HONAM BRANCH 34, Cheomdangangbyeon-ro 99beon-gil, Gwangsan-gu, Gwangju, Republic of Korea / T : 062-973-6626

 **VITSRO**
1522-6626
lmashop@naver.com
www.lmashop.com

NECK SKIN CARE
L.ma

THE VALUE OF LIGHT - VITSRO

INTRODUCTION OF VITSRO

LED PRO VITSRO is a company that creates beautiful values with light.

VITSRO is a LED-specializing company equipped with professional human resources and superior technology.

We are doing our best to develop products optimized for skincare based on our business philosophy of “creating beautiful values with light” and on LED technology.

We would like to present L.ma, which is a line of beauty devices for intensive care for the skin of the neck, in order to provide more customers with a beauty and happy life with LED light.

We promise to provide our customers with a more beautiful light through continuous development for higher values of LED.




LED WAVELENGTHS AND EFFECTS

The LED MASK is a product that use s phototherapy to treat damaged skin with light, according to the principle of promoting biochemical reactions in the skin.

According to research on the effect of a 660-nm wavelength on human skin cells, it improves skin elasticity by increasing the synthesis of collagen, which maintains skin elasticity, and by reducing the amount of MMP-1, which is an enzyme of decomposing collagen.



It uses a red wavelength that makes the skin healthy and elastic, a blue wavelength that relieves skin problems, and near-infrared rays that revitalize the skin cells by penetrating deep into the skin.



Near-infrared rays improve skin elasticity and revitalize skin cells by penetrating into the skin's subcutaneous tissue. Near-infrared rays, which were used by NASA (National Aeronautics and Space Administration) for phototherapy, are applied to be used by the beauty device. Invisible near-infrared rays are used together with the light of a red wavelength.



The light of a red wavelength synthesizes collagen and improves elasticity. It generates collagen by revitalizing skin cells after penetrating into the skin's dermal tissue and treats inflammation and wounds while making the skin appear youthful.



The light of a blue wavelength reduces skin problems and improves the skin. It reduces skin problems with its skin-soothing effect while preventing skin damage.

ADVANTAGES OF L.MA LED NECK MASK

ADVANTAGES OF L.MA LED NECK MASK

LED MASK L.MA DEDICATED TO THE SKIN OF THE NECK

-
-
-

Super light and practical

Embodies a light weight of 130g that does not feel heavy on the neck

Specialty

Reliable LED mask made by a Korean LED-specializing company

Availability

Premium LED mask at reasonable prices

Fantastic portability

Its light and small volume allows it to be used easily anytime, anywhere


Comfortable wear

It is easy to wear due to its back-side open shape. Its light weight and comfortable wear allows for light exercise and daily activity during use.

Sophisticated design


Anyone can use the device thanks to its refined design, regardless of age or gender.

Wavelengths suitable for the skin on the neck



The near-infrared rays, light of a red wavelength, and light of a blue wavelength are suitable to improve the skin on the neck. The LED at-home care product brightens the skin, reduces wrinkles, improves the skin density and elasticity, and soothes the skin with 108 LED lights suitable for the skin on the neck.


Convenient use



The action mode of L.ma can be selected with one button.
The device can be operated easily with just one power button.

L.MA CLINICAL TEST RESULTS


L.ma is a safe and effective product according to the results of a clinical test conducted by the clinical test agency.




The results of a clinical test on the neck wrinkles of 20 females ages 20s-50s demonstrated a 37% improvement rate after 4 weeks of product use.


Test Agency	Korea Testing & Research Institute, Healthcare Research Center, Skin Research Center	Test Subjects	20 females (ages 22-54)
Test Name	Clinical test on effects of VITSRO L.ma LED skincare device's reduction of neck wrinkles and improvement of skin tone	Test Method	Skin measurement before product use, skin measurement after 2 weeks and 4 weeks of product use, evaluation of abnormal reactions, survey evaluation

Product satisfaction survey



Improvement compared to the skin condition before and after 4 weeks of use





-
-
-

VITSRO works in a clean, reliable environment, and produces all of its products himself.

Based on our faith, we value all customers who have purchased VITSRO's products as if they were one of our family members, and we promise that we will faithfully carry out diverse guidance including.

We will work harder to learn more about the beautiful value of light, deliver it, and become a meaningful company.