

nelis

ENDO KEEPER



Endo keeper **CG series**

Endo Keeper CG series enables surgical instruments through the port (5-12mm) and hands to be put in the body and take out large specimens.



- **Simple**

Simple design



- **Detachable**

Can be opened during operation




- **Convertible**

Easy to convert laparoscopy to open surgery



ENDO KEEPER
CG-Type

ENDO KEEPER
Original-Type



Endo keeper CG series characteristics

Easy insert
with two step design



Can be opened
during operation
Easy to extract with
surgical specimen
during operation

- CG series

CO2 injection is possible
without the help of other instrument



| Endo Keeper Model List & Measurement

		Incision	Diameter	Length
Xsmall	Xsmall	20mm	41mm	130mm
Small	Small-S	20-40mm	60mm	185mm
	Small-S-C			265mm
	Small-S-CG			185mm
	Small-NS			265mm
	Small-NS-C			
	Small-NS-CG			
Medium	Medium-S	25-60mm	95mm	185mm
	Medium-S-C			265mm
	Medium-S-CG			185mm
	Medium-NS			265mm
	Medium-NS-C			
	Medium-NS-CG			
Large-N	Large-N	50-90mm	125mm	190mm
XXLarge	XXLarge	90-150mm	195mm	210mm

* S=string, NS=No string, C= Cap, G=Gas Valve



NELIS CO., LTD

1005, 201-dong, 397, Seokcheon-ro, Bucheon-si, Gyeong-gi-do, 14449, Korea

Tel : + 82 32 624 1697 **Tax :** + 82 32 624 1699

Email : nelismed@daum.net / nelismedical@naver.com

Website : nelis.co.kr



MEDITECH-INFRAMED CO., LTD

515, Dongnam-ro, Gangdong-gu, Seoul, Republic of Korea

Tel : + 82 2 476 8381 **Fax :** + 82 2 476 1640

Email : inframed@empas.com

website : inframed.co.kr