


REGENIN™ Filler

Want To Go Back?




NineBiopharm Co., Ltd.

New Generation Peptide Filler

REGENIN™

Regeneration from inside


- **Peptidomimetic: Beyond Peptide**



- **Immediate visible results
& Long term benefit**



- **Regeneration from inside**
 - increase collagen synthesis
 - decrease elastase activity



- About Us

Company at a glance

► NineBiopharm's ARisT platform


- NineBiopharm's proprietary ARisT (Advanced Remodeling in silico Technology) platform is a revolutionary hit finding the novel biomaterials
- NineBiopharm has new biomaterials that can be used for cosmetics, medical devices, and medicines
- NineBiopharm's new biomaterials can be used an effective treatment for various diseases
- In Korea, NineBiopharm is the most advanced biomaterials development company and moving fast toward commercialization and successful launches in markets across the world



• Differentiation of ReGENIN™ Filler

	J Filler	E Filler	ReGENIN™
pH	6.7~7.8	6.0~7.5	6.8~7.2
Residual Cross -linker (ppm)	≤2 (After sterilization)	≤2 (After sterilization)	≤2 (Before sterilization)
Injectability (N)	6~14	≤40	10~22
Osmolality (mmol/kg)	270~390	250~350	270~330
Endotoxin (EU/MI)	<0.5	<0.5	<0.5
HA contents	20mg/ml	20mg/ml	24mg/ml
Active Ingredient	X	X	Two-Peptidomimetic materials*
Price			

***Peptidomimetic:** Peptides or compounds with improved biological stability/activity, etc. through modification of the existing sequence and structure as a generic term for small substances designed to mimic the function of peptides




• Indication and Summary Table

REGENIN™					
Product Name	LIGHT	CLASSIC	ULTRA	VOLUME	LIPS
HA Concentration	24 mg/mL				
Visco-elasticity	Low	Mid	High	Hard	High
Packaging	1.1mL/Syringe				
Lidocaine	0.3%				
Injection Area	<ul style="list-style-type: none"> - Superficial wrinkles - Filling of all moderately pronounced wrinkles on the face - Eyelid Tear Troughs, Fine Lines, Upper Cheek, Hollow Cheek 	<ul style="list-style-type: none"> - Deep wrinkles - Harmonization of volume - Fore Head, Nose Bridge, Marionette Lines, Tear Troughs, Hollow Cheek 	<ul style="list-style-type: none"> - Very Deep wrinkles - Fore Head, Nose Bridge, Marionette Lines, Hollow Cheek, Nasolabial Folds 	<ul style="list-style-type: none"> - Creation of the volumizing & contouring - Jaw Lines, Chin, Mandibular Region, Tip of nose 	<ul style="list-style-type: none"> - Creation of the volumizing & contouring - Lips
Implant Duration Before touch up ※	12 months on average				
Shelf Life	24 months				
Storage	Store in 1-30°C. Protect from direct sunlight and freezing				



※ Average duration depends on several factors: Patients' age, skin type, lifestyle, severity of the wrinkles to be corrected, type of injection, the injected volume and the skill of a practitioner

• Applying Area





- In vivo data

- Roughness, Wrinkle Evaluation



- Collagen and Elastin level



• Product



REGENIN™ LIGHT

New Generation Peptide Filler Sodium Hyaluronate Gel with Lidocaine

1.1ml/Syringe

- Superficial wrinkles
- Filling of all moderately pronounced wrinkles on the face
- Eyelid Tear Troughs, Fine Lines, Upper Cheeks, Hollow Cheek



REGENIN™ CLASSIC

New Generation Peptide Filler Sodium Hyaluronate Gel with Lidocaine

1.1ml/Syringe

- Deep wrinkles
- Harmonization of volume
- Fore Head, Nose Bridge, Marionette Lines, Tear Troughs, Hollow Cheek



REGENIN™ ULTRA

New Generation Peptide Filler Sodium Hyaluronate Gel with Lidocaine

1.1ml/Syringe

- Very Deep wrinkles
- Fore Head, Nose Bridge, Marionette Lines, Cheeks, Hollow Cheek, Nasolabial Folds



REGENIN™ VOLUME

New Generation Peptide Filler Sodium Hyaluronate Gel with Lidocaine

1.1ml/Syringe

- Creation of the volumizing & contouring
- Jaw Lines, Chin, Mandibular Region, Tip of Nose



REGENIN™ LIPS

New Generation Peptide Filler Sodium Hyaluronate Gel with Lidocaine

1.1ml/Syringe

- Creation of the volumizing & contouring
- Lips

NineBiopharm Co., Ltd.

ReGENIN™ Filler



NineBiopharm Co., Ltd.

125, Osongsaengmyeong 2-ro, Osong-eup, Heungdeok-gu, Cheongju-si,
Chungcheongbuk-do, Republic of Korea

Tel: +82-43-903-9651 Fax: +82-43-903-9653

E-mail: nbpharm@ninebiopharm.com

Website: www.ninebiopharm.com