


Wireless Massager




UCARETRON

Who use it?



Workers



Housewives



Students



Drivers




Seniors

EASY CONTROL

The massager is connected to a smart phone by Bluetooth and controlled by app(Android/iOS)

As the massager is often put on the shoulder or back, wireless control is inevitable



STEREO

User can operate two massagers simultaneously

They can get a massage on different parts of their body at the same time just like hand-massage



Each massager is separately controlled by time-delay, stimulation level.



Music Beat Massage

The massager stimulates to the music beat.




40 songs are pre-saved in the app

- Intensive beat music : 20 songs
- Piano : 10 song
- Meditation music : 10 songs

The user can make a playlist with their own music




Massage Mode



Various Pad



Various Pad




HYDROGEL ON PAD

Hydrogel on the pad
maximizes effectiveness of
electrical stimulation.



If the adhesion is weak, it can be used
more after being washed in the water

Battery



PRODUCT COMPOSITION



SPECIFICATION

Model No	UM-200B	
Input	5V	
Current Consumption	20mA	
Battery	Li-ion Polymer, 3.7V / 70mAh	
Material	ABS / Polyurethane(Pad)	
Size(mm)	Beat Connected: 40x40 <ul style="list-style-type: none">• Small Pad : 146x57• Large Pad : 195x68• Asymmetric Pad : 177x70,• Wrist Pad : 175x35	
Weight	Beat Connected : 15.2g <ul style="list-style-type: none">• Small Pad : 14g• Large Pad : 27.6g• Asymmetric Pad : 22.6g• Wrist Pad : 14.2g	