

Korean Best Contact Lens Solution

UNIVERSAL Co., Ltd. The Best OEM Manufacture



TEL(office): +82-62-956-3553 (HP) +82-10-2036-6283. FAX:+82-62-956-3554



Free to work with Universal as No.1 Korean Contact lens Solution Manufacture! Universal can provide you with best price and best quality! Universial is factory itself, Every types of OEM is ok!

UNIVERSAL Co., Ltd			
CEO	Nam Byeong-Sook		
Address	682-2 Bia-dong, Gwangsan-gu, Gwangju, 506-300, Korea TEL(office): +82-62-956-3553, (HP): +82-10-2036-6283, Fax +82-62-956-3554 E-mail: securitad@hanmail.net		
Establishment	2013		
Manufacture	Contact lens Solution, Saline solution, over the counter(OTC)		
Business type	Factory, Manufacturer		



Friendly Company Delivering Value to Customers in Everyday Life with Product Development and Distribution Innovation!

Universal, as a professional OTC(over the counter) product manufacturer & distributor which carefully selects and provides a variety of high quality and lower price-products with customer's satisfaction as sharing vision with customers.

Universal always puts most important thing to transparent, respecting-people, natural- friendly management in practice.

Universal is customer - oriented and friendly business and always pursues to manufacture best –quality products and distribution innovation to give customer values.

From now on, we, **Universal** shall declare to always focus on social responsibility and obligation from the old-ethic thinking that we are one of the social member loved by our fellow customers ever.

April 24th, 2013 Nam Byeong Sook, CEO Universal

The Path of Universal has followed

2013 Feb. Founded a Universal as over the counter manufacture

Apr. Got a manufacturing certificatation(Over the counter) from KFDA

Apr. Launched Contact lens Solution "EYE Me, LENS Me" and Saline Solution "Beauty Linse"

Apr. Entered into a contract with Davich optical service chain Entered into a Contract with Orange optical service chain Entered into a Contract with Lens me optical service chain Contracted with sole distributor total 9 companies.

2015 May. Sold in domestic convenient stores such as NH HANARO Mart

Aug. Established Gyeongin, Yeongnam, Choongcheong Branchs
2017 Sep. Entered into a contract Vassen Vietnam(the biggest solution distributor in Vietnam)









INTO UNIVERSAL FACTORY - 100% Made in Korea















The most important thing is Customer satisfaction, what's more? Nothing.

Special Business Know-how

A lot of Business Experience on a variety of business fields

Prominent abilities to RD team to develop, manufacture the goods.

Manufacturer and Exporter of contact lens solution. Gained the acknowlegement by RD team We are one of the biggest Korean distributor in Korea,

To target oversea business partners

Exporting to Asia, Europe, North and South America


Ability to customize and Make OEM process

Cusomization and OEM of contact lens colution following the customer request.





Business Strategy







"Lens Me"



Brand Name Lens Me & Eye Me Manufacturer **Universal Corp** Classification OTC(Over the counter) Weight 40ml, 60ml, 80ml, 120ml 150ml, 350ml Main Ingredients Polyhexamethylene Biguanide 20% For Cleansing, Sterilizing, Rinsing, Usage Moisturizing, Removing protein of lens Store vacuum-sealed at room Storage temperature Valid Through 3 years from the date of manufacture





Universal Series Unversal Series-2 "Eye Me"



LENS Me, Eye Me: Multi-purpose solutions for waching, disinfecting and preserving at once!

- To make for your healthy eyes, To use to prevent the lens deformation!
- Powerful cleansing foreign substances such as protein, suplus fat, calsium which are easily stained with lens starting
 - To prevent the lens deformation and extend the lens durablity.
- Sterilizing "acanthamoeba cyst" bacteria
- Adding HPMC(Hydroxypropy Methylcellulose) for enhancing moisture on eyes.
 - To keep comfortable feeling and reduce the friction between eye and lens when you wear. (combining moisture and viscosity characteristic)
 - To keep protective moisture coating over the surface of Lens and pass the moisture constantly.
- Keeping proper acidity of tear level.
 - To Keep comfortable feeling as if you wear for the first time by maintaining water content and form, dimension, Sagittal value. (Stable PH level)
 - To relieve tingly and chilled feeling absolutely and suppress fatigue of eyes.



Main Components – Multi purpose solution



Lense case is optional!



MAIN COMPONENT	HOW TO ACT	
EDTA	To prevent from Lens deformation	
РНМВ	To sterilize Acanthamoeba bacteria	
Ploxamer 407	To remove protein of lens	
НРМС	Artificial tear – Moisturizing and prevention from dry-feeling on eyes	
Boracic acid	To keep PH 6.0 ~ 8.0	
Sodium chloride	To control osmotic pressure	



Main Product line



Multi-Purpose lens solution

- MPS(Multi-purpose solution) for disinfecting acanth amoeba by using main harmless element PHMB(Poly hexa methylene Biguanide)
- No need to other process of protein removal / cleaner / salt solution(heat)disinfection

Solution(40, 60, 80, 120, 150, 350ml)



Drop O2(13ml)



- Intensifying moisture by using HPMC(Hydroxy Propy Methyl Cellulose)
- Keeping comfortable wearing by maintaining humidity on lens surface

No worry for Yellow sand and traveling

- Combining all benefits : Washuing + Moisturizing + Safe Preserving
- Protable small-sized products in yellow-sand or travel.



Beauty lense(500ml)



Make your own brand!!







Around 80~90 Optical brands over the worlds OEM from Universal!!










Make your own design!!

Any languages and contents are supported for your insert and box design
















Fourfold & Deproteinize system



Step 1TRIKLENS of fourfold deprotenizing system moves to remove protein deposits on the surface of lens



Step 3 A surfactant, Poloxamer 407 separates protein deposits and push up from the surface of lens.



Step 2 EDTA and Phosphate Buffer break off calcium connectors supporting protein deposits.



Step 4Protein deposits separated and removed from lens by triple deproteinizing system.



How to use the contact lens solution?



How to Use (Step)

- 1. Squirt contact lense solution drop on the contacts in the palm of your hand (You can use lense case or clean surface of table)
- 2. Take your finger for a few seconds and slowl rub your contacts in circular motions. (Rub for 10 seconds)
 After that, the high resolution of contact and contact other hand rub
- 3. Leave it as they are for more than 5 minutes (If possible, We suggest you to use Lense case)
- 4. Once that is done, put it in your contacts. Contact your submersed in liquid solutions that ensure that enough.

Efficiency

- 1. 4-hours Care(Rinsing without rubbing)
- 2. Place cleaned lenses in lens case and fill with this Solution.
- 3. Soak for at least four hours with lids close
- 4. Expiry after 3 Years







Compare with other products

	Universal	Bausch&Lomb ReNu	Alcon OptiFree
Main Ingredients	20% polyhexamethylene biguanide hydrochloride PHMB 0.1mg (EDTA 0.025%Pluronic F127) Polyquatrinum-1 3.3mg EDTA Na 50mg	PHMB 0.05mg (EDTA 0.12% Tetronic)	Tetradecanamide 0.5mg Polyquatrinum-1 3.3mg EDTA Na 50mg
Regimen	5 minutes 4 hours (12 hours Moist)	4 hours	No rub – 6 hours
Storing	30days	-	30days
Lens Case	NO	No	No




Certificated by Korean Government









- Ingredient Table -

Thank you