

# MaruAra Color Lens

Our Circle Lens maintains healthy view and manage clean eyesight of the people of the whole world and expresses naturalness and beauty of the eyes. Our lenses are easy to wear, comfortable and available to anyone without age restrictions. It also protects UV rays from at least 98% to up to 99%, reducing eye fatigue.



### Silicone Hydrogel Lens


#### Eye View Silicone Hydrogel Lens offer the followings

Silicone hydrogel lens adapt the material to have advantages such as high oxygen permeability of hard lens, wearing sensation of soft lens, and the percentage of water content by unifying the advantages of hard lens and soft lens as one.

- Oxygen permeability(DK) : 27×10<sup>11</sup>(cm2/sec)(mlO2/ml×mmHg) Water Content : 52% - D I A : 14.2mm B . C : 8.6mm



#### Daily wearable Soft Contact Lens Uniform quality from cast molding method of production THE LENS lens made of hydrophilic material is soft and flexible. It boasts of




#### Clear refracting power and good wetting force


Our lens boasts of clear refractive power due to its high water absorption function and excellent optical properties. In addition, it has good tears wetting force with moisturing power.

#### FDA approved cosmetic pigments


THE LENS lens provides superior pigment stability by inserting cosmetic pigment into the material of lens which was approved by FDA, satisfies the customer's requirements with a variety of designs and pigment tone.




THE LENS's has excellent shape stability so that the beginners of lens wearers can be easily adapted to the lens, and provides enhanced stability by minimizing irritation to the corneal surface and the eyelids due to the thin and finely treatment of edges.




#### **UV** Test




# MaruAra Lens Edge















N Brand Edge


жыс - 461.26 Цара - 105





A Brand Edge





B Brand Edge

## MaruAra Circle Lens Color

2019 Design







Gray Color



Choco Color


Brown Color

# MaruAra Circle Lens Color





We always Respect and cherish the wellness and beauty of you. To satisfy your high leveled demand. We released new products. Enoy The Lens Contact Lens.









We always Respect and cherish the wellness and beauty of you. To satisfy your high leveled demand. We released new products. Enoy The Lens Contact Lens.






demand. We released new products. Enoy The Lens Contact Lens.


GNE













We always Respect and cherish the wellness and beauty of you. To satisfy your high leveled demand. We released new products. Enoy The Lens Contact Lens.












BLUE HSD-44









B.C	8.2mm ~ 8.8mm
DIA	13.5mm ~ 15.5mm
Power	±0~±25
Materials	2HEMA+a
Water Contents	38%~55%
Color	Gray,Brown,Green,Blue,Violet,Aqua,Hazel



# Package





Rebecca	Rebecca
CONTENTS the Set Constitution A contents of the Set Constitution A before Set Constitution and Set The Net Constitution and Set The Net Constitution and Set Set Set Set Constitution and Set Set Set Set Set Set Set Set Set Set Set Set Set Set	CONTROLS ON BUTCH IN A
0 1	2







The manual is located inside the box.





Address: #401, 311-8, Tanjung-ro, Ilsandong-gu, Goyang-si, Gyeonggi,do, Korea Zip Code: 10336 Web-Site: http://maruara.quv.kr E-Mail: Sky J. Lee <u>maruara0620@gmail.com</u> Tel: 031-913-8628, 031-977-9828 Fax: 031-977-9829