



Skinable is a combination of
Skin and -able,
meaning everyone can have beautiful skin.
We are making skin care device
for natural skin improvement
aiming at low irritation skin care.









Galvanic Massager True Solution LED Massager Face Cooler

Galvanic Massager





Sebum

Removal





Nutrient

Absorption



Pore Elasticity

<u></u>

Triple-care(peeling, nourishment, elasticity) with a single device

24K pure gold blade safe for skin

30,000-time vibrations per second

Effects of Galvanic

Using the feature of identical electrodes pushing away each other, it helps the absorption of nourishing ingredients in cosmetic products that are not easily absorbed by patting with hands by transforming such ingredients into cations or anions.



Anion Penetration of nutrients



True Solution LED Massager







Galvanic

Cation · Anion



Microwave

Vibration



LED

Therapy

Ultrasound Massage

cellular tissues

Circulation improvement and wrinkle care through increase in nour temperature of

Polar
ionization vit
pulls away dirt
and helps c
nourishing w
ingredients
permeate
deep into skin

200-time Cell vibrations per second in diminishes wrinkles and steri improves elasticity im

Cell regeneration, acne and inflammation reduction, sterilization effect, skin tone improvement, melanin suppression

Effects of LED Lights



Rays penetrate up to dermal tissue of the skin Helps skin elasticity and regeneration by helping collagen production Rays used to treat inflammation, wounds and hair loss



Rash, sensitive skin and oily skin improvement effects

Effective in improving skin troubles or acne symptoms by eliminating bacteria that cause inflammation



Effective in diminishing fine lines
Helps skin metabolism for vitality
Skin calming, skin tone improvement and suppression of excessive melanin

Face Cooler











Reduction of Reduction of Heated Puffiness Skin Temperature

Skin Calming

Pore Elasticity

Differentiated cute design that can be used as an accessory holder

Replaceable use with separable coolant

Functionally designed for the area that comes in contact with skin

Long-lasting cooling sensation with harmless refrigerant

Recommended for the Following

To soothe skin after being exposed to heat

To tighten pore after basic skin care steps

To de-puff face and legs

To relieve of various pains such as tooth pain, mastitis To alleviate swelling after cosmetic surgery or laser treatment

To relieve swelling caused by poor blood circulation

For relieving skin irritated after depilation When face becomes swollen after eating at night

te When in need of skin relief ter before t putting on makeup

As the finishing stage after use of Skinable massager